



innovation for health & wellness

“trace & catch”

Instructions for Use

Certolizumab pegol (Cimzia®) ELISA

SHIKARI® Q-CERT

Enzyme immunoassay for the quantitative determination of certolizumab pegol (Cimzia®) in serum and plasma

REF TR-CERTv1



12 x 8



2-8 C

Revision # 1.1 August 2018



Matriks Biotek® Laboratories
www.matriksbiotek.com



Contents	Page
Intended Use	3
Summary and Explanation.....	3
Test Principle	4
Warnings and Precautions.....	4
Storage and Stability.....	5
Specimen Collection and Storage	5
Materials Supplied.....	6
Materials Required but not Supplied	6
Procedure Notes.....	7
Pre-Test Setup Instructions.....	8
Test Procedure.....	9
Quality Control	10
Calculation & Interpretation of Results	10
Assay Characteristics	12
Automation	12
References.....	13
Semi-Log Graph Paper	14
Semi-Log Graph Paper	15

	SHIKARI [®] Q-CERT
	Free certolizumab pegol (Cimzia [®]) quantitative analyses
Required Volume (µl)	10
Total Time (min)	70
Sample	Serum, plasma
Sample Number	96
Detection Limit (ng/mL)	12
Spike Recovery (%)	Between 85-115
Shelf Life (year)	1

Intended Use

Enzyme immunoassay for the quantitative determination of **free certolizumab pegol** (Cimzia®) in serum and plasma. *Matriks Biotek® Certolizumab pegol ELISA* has been especially developed for the quantitative analysis of free certolizumab pegol in serum and plasma samples at high specificity.

Summary and Explanation

Certolizumab pegol is an anti-TNF (anti-tumour necrosis factor) drug. It blocks the action of TNF protein and so reduces inflammation in conditions including rheumatoid arthritis, psoriatic arthritis and ankylosing spondylitis.

TNF α is a key pro-inflammatory cytokine with a central role in inflammatory processes. Biological activity associated to TNF α include the upregulation of cellular adhesion molecules and chemokines, upregulation of major histocompatibility complex (MHC) class I and class II molecules, and direct leukocyte activation. TNF α stimulates the production of downstream inflammatory mediators, including interleukin-1, prostaglandins, platelet activating factor, and nitric oxide. After treatment with certolizumab pegol, patients with Crohn's disease demonstrated a decrease in the levels of C-reactive protein (CRP).

Certolizumab pegol binds to free and membrane-bound human TNF α with a KD of 90pM and neutralizes its activity. Extent of neutralization is also dose-dependent. It also inhibited the release of lipopolysaccharide-induced IL-1 β from monocytes. TNF α is a key pro-inflammatory cytokine with a central role in inflammatory processes in which elevated levels have been observed in patients with RA and Crohn's. Certolizumab pegol selectively neutralizes TNF α . It does not bind to TNF- β . As certolizumab is only a Fab' fragment and thus missing the Fc region, it does not fix complement or cause antibody-dependent cell-mediated cytotoxicity. Furthermore, apoptosis of monocytes or lymphocytes, or neutrophil degranulation have not been observed in vitro.

There is a linear relationship between dose administered and Cmax and AUC. A mean Cmax of approximately 43 to 49 mcg/mL occurred at Week 5 during the initial loading dose period using the recommended dose regimen for the treatment of patients with rheumatoid arthritis (400 mg sc at Weeks 0, 2 and 4 followed by 200 mg every other week). Tmax, SubQ dose = 54 - 171 hours; Bioavailability, SubQ dose = 80% (range of 76% - 88%)

Test Principle

Solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle. Standards and samples (serum or plasma) are incubated in the microtitre plate coated with the reactant for certolizumab pegol (Cimzia®). After incubation, the wells are washed. A horse radish peroxidase (HRP) conjugated probe is added and binds to certolizumab pegol captured by the reactant on the surface of the wells. Following incubation wells are washed and the bound enzymatic activity is detected by addition of chromogen-substrate. The colour developed is proportional to the amount of certolizumab pegol in the sample or standard. Results of samples can be determined directly using the standard curve.

Warnings and Precautions

1. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood. For further information (clinical background, test performance, automation protocols, alternative applications, literature, etc.) please refer to the local distributor.
3. In case of severe damage of the kit package please contact Matriks Biotek® or your supplier in written form, latest one week after receiving the kit. Do not use damaged components in test runs, but keep safe for complaint related issues.
4. Obey lot number and expiry date. Do not mix reagents of different lots. Do not use expired reagents.
5. Follow good laboratory practice and safety guidelines. Wear lab coats, disposable latex gloves and protective glasses where necessary.
6. Reagents of this kit containing hazardous material may cause eye and skin irritations. See MATERIALS SUPPLIED and labels for details.
7. Chemicals and prepared or used reagents have to be treated as hazardous waste according the national biohazard safety guidelines or regulations.
8. Avoid contact with Stop solution. It may cause skin irritations and burns.
9. All reagents of this kit containing human serum or plasma (i.e. standards) have been tested and were found negative for HIV I/II, HBsAg and HCV. However, a presence of these or other infectious agents cannot be excluded absolutely and therefore reagents should be treated as potential biohazards in use and for disposal.
10. Some reagents contain sodium azide (NaN_3) as preservatives. In case of contact with eyes or skin, flush immediately with water. NaN_3 may react with lead and copper plumbing to form explosive metal azides. When disposing reagents, flush with large volume of water to avoid azide build-up.

Storage and Stability

The kit is shipped at ambient temperature and should be stored at 2-8°C. Keep away from heat or direct sun light. The strips of microtiter plate is stable up to the expiry date of the kit in the broken, but tightly closed bag when stored at 2-8°C.

Specimen Collection and Storage

Serum, Plasma (EDTA, Heparin)*

The usual precautions for venipuncture should be observed. It is important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Do not use grossly hemolytic, icteric or grossly lipemic specimens. Samples appearing turbid should be centrifuged before testing to remove any particulate material.

Storage:	2-8°C	-20°C	Keep away from heat or direct sun light Avoid repeated freeze-thaw cycles
Stability:	2 d	6 mon	

*. Certolizumab pegol (Cimzia®) infusion camouflages/masks the presence of antibody to certolizumab pegol (ATC) in serum/plasma samples. Therefore, blood sampling time is critical for detection of ATC. The Matriks Biotek® Laboratories suggests to obtain blood sample just before the infusion of certolizumab pegol (Cimzia®) or at least 2 weeks after the infusion of certolizumab pegol (Cimzia®)

Materials Supplied

1 x 12 x 8	MTP	Microtiter Plate Break apart strips. Microtiter plate with 12 rows each of 8 wells coated with reactant.
8 x 0.5 mL	STND A-F HIGH CNTRL LOW CNTRL	Certolizumab pegol Standards A-F, High Level Control, Low Level Control 1000; 333; 111; 37;12; 0 nanogram/mL Ready to use. Used for construction of the standard curve. Contains certolizumab pegol (Cimzia®), human serum, stabilizer and <0.1% NaN ₃ .
1 x 50 mL	ASSAY BUF	Assay Buffer Blue colored. Ready to use. Contains proteins and <0.1% NaN ₃ .
1 x 12 mL	HRP CONJ	Horse radish peroxidase-Conjugated Probe Red colored. Ready to use. Contains HRP-probe, stabilizer and preservatives.
1 x 12 mL	TMB SUBS	TMB Substrate Solution Ready to use. Contains TMB
1 x 12 mL	TMB STOP	TMB Stop Solution Ready to use. 1N HCl.
1 x 50 mL	WASHBUF CONC	Wash Buffer, Concentrate (20x) Contains Buffer with Tween 20.
2 x 1	FOIL	Adhesive Foil For covering of Microtiter Plate during incubation.
2x 1	SLGP	Semi-Log Graph Paper For constructing standard curve and calculation of results.

Materials Required but not Supplied

1. Micropipettes (< 3% CV) and tips to deliver 5-1000µL.
2. Calibrated measures.
3. Tubes (1 mL) for sample dilution.
4. Wash bottle, automated or semi-automated microtiter plate washing system.
5. Microtiter plate reader capable of reading absorbance at 450/650 nm.
6. Bidistilled or deionised water, paper towels, pipette tips and timer.

Procedure Notes

1. Any improper handling of samples or modification of the test procedure may influence the results. The indicated pipetting volumes, incubation times, temperatures and pretreatment steps have to be performed strictly according to the instructions. Use calibrated pipettes and devices only.
2. Once the test has been started, all steps should be completed without interruption. Make sure that required reagents, materials and devices are prepared ready at the appropriate time. Allow all reagents and specimens to reach room temperature (18-25 °C) and gently swirl each vial of liquid reagent and sample before use. Mix reagents without foaming.
3. Avoid contamination of reagents, pipettes and wells/tubes. Use new disposable plastic pipette tips for each reagent, standard or specimen. Do not interchange caps. Always cap not used vials. Do not reuse wells/tubes or reagents.
4. Use a pipetting scheme to verify an appropriate plate layout.
5. Incubation time affects results. All wells should be handled in the same order and time sequences. It is recommended to use an 8-channel Micropipettor for pipetting of solutions in all wells.
6. Microplate washing is important. Improperly washed wells will give erroneous results. It is recommended to use a multichannel pipette or an automatic microplate washing system. Do not allow the wells to dry between incubations. Do not scratch coated wells during rinsing and aspiration. Rinse and fill all reagents with care. While rinsing, check that all wells are filled precisely with Wash Buffer, and that there are no residues in the wells.
7. Humidity affects the coated wells/tubes. Do not open the pouch until it reaches room temperature. Unused wells/tubes should be returned immediately to the re-sealed pouch including the desiccant.

Pre-Test Setup Instructions

1. Preparation of Components

Dilute/ dissolve	Component	with	Diluent	Relation	Remarks	Storage	Stability
10 mL	Wash Buffer*	Up to 200 mL	bidist. Water	1:20	Warm up at 37°C to dissolve crystals. Mix vigorously.	2-8 °C	2 w

*. Prepare Wash Buffer before starting assay procedure.

2. Dilution of Samples

Sample	To be diluted	With	Remarks
Serum/ Plasma	1:50	Assay Buffer	For dilution at 1:50; 10 µl Sample + 490 µl Assay Buffer

Patient samples with a concentration of certolizumab pegol above the measuring range are to be rated as > "Highest Standard (Standard A)". The result must not be extrapolated. The patient sample in question should be further diluted with Assay Buffer and retested.

Test Procedure

1	Pipette 50µl of Assay Buffer non-exceptionally into each of the wells to be used.
2	<p>Pipette 50 µL of each ready-to use Standards, High Level Control, Low Level Control and Diluted Samples into the respective wells of microtiter plate.</p> <p>Wells</p> <p>A1: Standard A B1: Standard B C1: Standard C D1: Standard D E1: Standard E F1: Standard F G1: High Level Control H1: Low Level Control A2 and on:Sample (Serum / Plasma)</p>
3	Cover the plate with adhesive foil. Incubate 30 min at room temperature (18-25°C).
4	Remove adhesive foil. Discard incubation solution. Wash plate 3 times each with 300 µL of diluted Wash Buffer. Remove excess solution by tapping the inverted plate on a paper towel.
5	Pipette 100 µL of ready-to use HRP-Conjugated Probe into each well.
6	Cover the plate with adhesive foil. Incubate 30 min at room temperature (18-25°C).
7	Remove adhesive foil. Discard incubation solution. Wash plate 3 times each with 300 µL of diluted Wash Buffer. Remove excess solution by tapping the inverted plate on a paper towel.
8	Pipette 100 µL of TMB Substrate Solution into each well.
9	Incubate 10 min (without adhesive foil.) at room temperature (18-25°C) in the dark.
10	Stop the substrate reaction by adding 100 µL of Stop Solution into each well. Briefly mix contents by gently shaking the plate. Color changes from blue to yellow.
11	Measure optical density with a photometer at 450/650 nm within 30 min after pipetting of the Stop Solution.

Quality Control

The test results are only valid only if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable standards/laws. All standards must be found within the acceptable ranges as stated above and/or label. If the criteria are not met, the run is not valid and should be repeated. In case of any deviation the following technical issues should be proven: Expiration dates of (prepared) reagents, storage conditions, pipettes, devices, incubation conditions and washing methods.

Calculation & Interpretation of Results


1. Using the standards (1000; 333; 111; 37; 12; 0 ng/mL) disregarding zero standard, construct a standard curve by plotting the OD_{450/650} nm for each of 4 standards on the vertical (Y-axis) axis versus the corresponding certolizumab pegol concentration on the horizontal (X-axis) axis, thus creating a standard curve by 4 points obtained.
2. The concentration of the samples can be read directly from this standard curve. Using the absorbance value for each sample, determine the corresponding concentration of certolizumab pegol from the standard curve. Find the absorbance value on the Y-axis and extend a horizontal line to the curve. At the point of intersection, extend a vertical line to the X-axis and read the certolizumab pegol concentration for the unknown sample.
3. If computer data regression is going to be used, we recommend primarily "4 Parameter Logistic (4PL)" or secondly the "point-to-point calculation".
4. To obtain the exact values of the samples, the concentration determined from the standard curve must be multiplied by the dilution factor (50x). Any sample reading greater than the highest standard should be further diluted appropriately with Assay Buffer and retested. Therefore, if the pre-diluted samples have been further diluted, the concentration determined from the standard curve must be multiplied by the further dilution factor.

E.g.; If the pre-diluted sample further diluted in a ratio of 1:5 then results should be multiplied by 250.

5. Automated method: Computer programs can also generally give a good fit.
6. For the OD values of **High Level** and **Low Level** controls, please refer to **Quality Control Certificate** (QCC) provided by each kit

Typical Calibration Curve

(Example. Do not use for calculation!)



Standard	Concentration (ng/mL)	Mean OD450/650nm
A	1000	2,701
B	333	1,977
C	111	1,006
D	37	0,515
E	12	0,235
F	0	0,171

Assay Characteristics

1. **Specificity:** Except for certolizumab pegol, there is no cross reaction with other therapeutic antibodies and native serum immunoglobins.
2. **Sensitivity:** The lowest detectable level that can be distinguished from the zero standard is 12 ng/mL.
3. **Precision Of Kit:**
Intra-assay CV: <15% for certolizumab pegol range 12-1000 ng/mL.
Inter-assay CV: <15% for certolizumab pegol range 12-1000 ng/mL.
4. **Recovery:** Recovery rate was found to be between 85-115% with normal human serum samples with known concentrations.

Automation

Experiments have shown that the **Matriks Biotek**[®] SHIKARI[®] Certolizumab pegol ELISA is also suitable to run on an automated ELISA processor.

References

1. CIMZIA® [prescribing information], Smyrna, GA: UCB, Inc.; 2018.
2. Keystone E, van der Heijde D, Mason D Jr, et al. Certolizumab pegol plus methotrexate is significantly more effective than placebo plus methotrexate in active rheumatoid arthritis: findings of a fifty-two-week, phase III, multicenter, randomized, double-blind, placebo-controlled, parallel-group study. *Arthritis Rheum.* 2008;58:3319-3329.
3. Schwartzman S, Morgan GJ Jr. Does route of administration affect the outcome of TNF antagonist therapy? *Arthritis Res Ther.* 2004;6(Suppl 2):S19-S23.
4. Sheikhzadeh A, Yoon J, Formosa D, et al. The effect of a new syringe design on the ability of rheumatoid arthritis patients to inject a biological medication. *Appl Ergon.* 2012;43:368-375.
5. Data on file. UCB, Inc.; Smyrna, GA
6. Veronese FM, Mero A. The impact of PEGylation on biologic therapies. *Biodrugs.* 2008;22:315-329.
7. Weir N, Athwal D, Brown D, et al. A new generation of high-affinity humanized PEGylated Fab' fragment anti-tumor necrosis factor-alpha monoclonal antibodies. *Therapy.* 2006;3:535-545.
8. Chapman AP. PEGylated antibodies and antibody fragments for improved therapy: a review. *Adv Drug Deliv Rev.* 2002;54:531-545.
9. Harris JM, Chess RB. Effect of pegylation on pharmaceuticals. *Nat Rev Drug Discov.* 2003;2:214-221.
10. Mease PJ, Fleischmann R, Deodhar AA, et al. Effect of certolizumab pegol on signs and symptoms in patients with psoriatic arthritis: 24-week results of a phase 3 double-blind randomised placebo-controlled study (RAPID-PsA). *Ann Rheum Dis.* 2014;73:48-55.
11. Landewe R, Braun J, Deodhar A, et al. Efficacy of certolizumab pegol on signs and symptoms of axial spondyloarthritis including ankylosing spondylitis: 24-week results of a double-blind randomised placebo-controlled phase 3 study. *Ann Rheum Dis.* 2014;73(1):39-47.
12. Weinblatt ME, Fleischmann R, Huizinga TW, et al. Efficacy and safety of certolizumab pegol in a broad population of patients with active rheumatoid arthritis: results from the REALISTIC phase IIIb study. *Rheumatology.* 2012;51:2204-2214.
13. Gladman D, Fleischmann R, Coteur G, Woltering F, Mease PJ. Effect of certolizumab pegol on multiple facets of psoriatic arthritis as reported by patients: 24-week patient-reported outcome results of a phase III, multicenter study. *Arthritis Care Res.* 2014;66(7):1085-1092.
14. Schreiber S, Khaliq-Kareemi M, Lawrance IC, et al. Maintenance therapy with certolizumab pegol for Crohn's disease. *N Engl J Med.* 2007;357:239-250.
15. Schreiber S, Colombel JF, Bloomfield R, et al. Increased response and remission rates in short-duration Crohn's disease with subcutaneous certolizumab pegol: an analysis of PRECISE 2 randomized maintenance trial data. *Am J Gastroenterol.* 2010;105(7):1574-1582.
16. Feagan BG, Coteur G, Tan S, Keininger DL, Schreiber S. Clinically meaningful improvement in health-related quality of life in a randomized controlled trial of certolizumab pegol maintenance therapy for Crohn's disease. *Am J Gastroenterol.* 2009;104(8):1976-1983.
17. Pallis AG, Mouzas IA, Vlachonikolis IG. The inflammatory bowel disease questionnaire: a review of its national validation studies. *Inflamm Bowel Dis.* 2004;10(3):261-269.

